



## DNA Safe Stain

For Research Use Only

**Cat. No.:** EP5082

**Storage:** Store at 2-8 °C (protected from light). **Spin down before use.**

**Quantity:** 1ml

**Shipment:** Room temperature

### **Description:**

DNA safe Stain is a new nucleic acid stain which can be used as a safer alternative to the traditional ethidium bromide stain for detecting nucleic acid in agarose gels. It is as sensitive as Ethidium bromide and can be used exactly the same way in agarose gel electrophoresis.

DNA safe Stain emits green fluorescence when bound to DNA or RNA. It has two secondary fluorescence excitation peaks (~300nm; ~400nm) and one strong excitation peak centered around 500nm. The fluorescence emission is centered at ~540 nm. Thus, DNA safe Stain is compatible with a wide variety of gel reading instruments.

DNA safe Stain can be used for precast agarose gels and when better sensitivity is needed - poststaining is recommended.

**Note:** Warm DNA Safe Stain to room temperature before use.

### **Applications:**

Non-carcinogenic alternative to Ethidium bromide.

### **Safety:**

DNA safe Stain is non-carcinogenic and according to the Ames test it causes significantly fewer mutations than Ethidium bromide.

### **Protocol:**

#### **Precasting:**

Prepare 100ml of agarose gel solution (concentration from 0.8-3.0%) and heat until the solution is completely clear and no small floating particles are visible.

Add **1-1.5µl** of DNA safe Stain to the gel solution and mix it gently.

Cool the gel to 50-60 °C and cast the gel, into the gel tray.

When the gel is solid, load the samples and perform electrophoresis.

Detect the bands under UV illuminator.

#### **Poststaining:**

The DNA safe Stain poststaining solution may be used 2-3 times. Staining solution to be reused should be preferably stored at room temperature in the dark.

For <0.5 cm thick agarose gel, 5-15µl of the stain should be used per 100ml of buffer.

Optimal staining time (5 - 60 minutes) and the amount of the stain may depend on the thickness of the gel and the percentage of agarose.

**Notes:** 1 ml of DNA safe Stain is sufficient for 50-65 L of agarose gel. The thickness of gel should < **0.5cm**. Repeated melting of gels containing DNA safe Stain may result in low sensitivity. DNA safe Stain is non-carcinogenic but may irritate skin and eyes. Please wear gloves while handling.

**SinaClon BioScience**

[www.sinaclon.com](http://www.sinaclon.com)

Tel: +98(0)21 4463 0050-51

+98(0)21 4466 5156

Customer Service: +98(0)9023120059

+98(0)21 4466 5156

Place your order: [order@sinaclon.com](mailto:order@sinaclon.com)