

Prestained Protein Ladder,(10-170KDa)

Cat. No. SL7011 **Size:** 500 µl

Cat. No. SL7012 **Size:** 100 µl

- 3 µl or 5 µl per loading for clear visualization during electrophoresis on 15-well or 10-well mini-gel, respectively
- 2-3 µl per well for general Western transferring.
- Apply more for thicker (>1.5 mm) or larger gel.

Description:

The Prestained Protein Ladder is a three-color protein standard with 10 pre-stained proteins covering a wide range molecular weights for 10 to 180 kDa. Proteins are covalently coupled with a blue chromophore except for two reference bands (one green and one red band at 25 kDa and 75 kDa respectively) when separated on SDS-PAGE (Tris-Glycine buffer). The Prestained Protein Ladder is designed for monitoring protein separated during SDS polyacrylamide gel electrophoresis, verification of Western transfer efficiency on membranes (PVDF, nylon or nitrocellulose) and for approximate sizing of proteins. The ladder is supplied in gel loading buffer and is ready to use. Do not heat, dilute, add reducing agent before loading.

Contents:

Approximately 0.2~0.4 mg/ml of each protein in buffer (20 mM Trisphosphate pH 7.5 at 25°C), 2% SDS, 1 mM 2-Mercaptoethanol, 3.6 M Urea, and 15% (v/v) Glycerol).

Quality Control:

5 µl of Prestained Protein Ladder resolves 10 bands in 4-20% SDS-PAGE (Tris-Glycine buffer) and after Western blotting to PVDF membrane.

Sale Office : No. 56, Azimi St., Nafisi Ave., Phase 1, Ekbatan, Tehran, Iran

Tel : 02144666204 02144631124


Storage:

Stable for up to 2 weeks at 25°C.

Stable for up to 3 months at 4°C. For long term storage, store at -20°C.

Guide for Molecular Weight Estimation (kDa):

Migration patterns of Prestained Protein Ladder in different electrophoresis conditions are listed below:



Note. The apparent molecular weight of each protein (kDa) has been determined by calibration against an unstained protein ladder in each electrophoresis condition. * supplement data should be considered for more accurate adjustment.

All products are for research use only. Caution: Not intended for human or animal diagnostic or therapeutic uses.