


RNX- Plus

(Solution for total RNA isolation)

REF

EX6101

Quantity: 25 ml



(For long storage, store at -20°C)

RUO



Components

| Contents | Amounts |
|-----------|---------|
| RNX- Plus | 25 ml |

Description

RNX - Plus is a Guanidine/phenol solution for total RNA isolation from homogenized sample. Through the action of Guanidine salt in RNA isolation procedure, simultaneously DNA and protein are precipitated in phenol phase. Aqueous phase contains all types of the genomic RNA with high quality.

RNX - Plus solution designed to isolate total RNA from different amounts of biological material.

The obtained RNA is ready for use in all downstream applications like: RT-PCR, cDNA synthesis, Northern, dot, and slot blot analyses, Primer extension, Poly A+ RNA selection and etc.

Starting Materials

| Row | Materials | Amounts |
|-----|-----------------------------|---|
| 1 | Cell Culture | Up to 1×10^7 cells, depending on the cell line |
| 2 | Bacterial cell | Up to 1×10^7 cells |
| 3 | Tissue | 30-50 mg |
| 4 | Yeast cell | 5×10^7 cells |
| 5 | Plant and Filamentous Fungi | Up to 100 mg |
| 6 | Liquid materials like serum | 100 μ l |
| 7 | Whole blood | At least 500 μ l (up to 2ml) |

(RNX-008-00/01) (1)

*Prior to RNA extraction from whole blood, RBCs must be removed by RBC Lysis Buffer and white blood cells (WBC pellet) must use as starting material.

*Use fresh blood sample for RNA blood extraction.

Required materials and tools that are not supplied with the kit:

| Row | Materials and Tools | SinaClon Catalog |
|-----|---------------------|------------------|
| 1 | Chloroform | - |
| 2 | Isopropanol | - |
| 3 | Ethanol 75% | - |
| 4 | DEPC treated water | CH8141 |
| 5 | RBC Lysis Buffer | EX6122 |
| 6 | DNase I | MO5401 |
| 7 | DEPC | CH8131 |

Protocol for RNA isolation

- 1) Add 1ml ice cold RNX™ – PLUS solution to 2ml tube containing homogenized sample.
- 2) Vortex 5-10 sec and incubate at room temperature for 5min.
- 3) Add 200 μ l of Chloroform.
- 4) Mix well for 15 sec by shaking (Do not vortex).
- 5) Incubate on ice or 4 °C for 5 min.
- 6) Centrifuge at 12000 rpm at 4 °C for 15 min.
- 7) Transfer the Aqueous phase to new RNase-free 1.5ml tube, (do not disturb the mid phase) and add equal volume of Isopropanol.
- 8) Gently mix and incubate on ice for 15 min.
- 9) Centrifuge the mixture at 12000 rpm at 4 °C for 15 min.
- 10) Discard the supernatant and add 1 ml of 75% Ethanol, shortly vortex to dislodge the pellet and then centrifuge at 4 °C for 8 min at 7500 rpm.
- 11) Discard the supernatant and let the pellet to dry at room temperature for few minutes (do not let dry completely, it will decrease the solubility of the pellet).
- 12) Dissolve pellet in 50 μ l of DEPC treated water. To help dissolving, place the tube in 55-60 °C water bath for 10 min.

(2)

There are different methods

There are two distinct and essential steps for RNA isolation by RNX - Plus solutions: Disruption and homogenization of sample.

Insufficient disruption and homogenization significantly will reduce RNA yield.

| Sample | Disruption | Homogenization |
|-----------------------------|--|--------------------|
| Cells | RNX solution | Vortex |
| Tissue | By mortar and pestle liquid Nitrogen | Syringe and Needle |
| Yeast cell | lyticase | Vortex |
| Plant and Filamentous Fungi | By mortar and pestle in liquid Nitrogen | Shredder |
| Liquid materials | RNX solution | Vortex |
| WBC pellet | RNX solution | Vortex |

Precautions





RNX - Plus contains an irritant (Guanidine thiocyanate) and poison (phenol). Handle with gloves and do not get in eyes, skin or clothing. Avoid breathing vapor.

In case of contact Immediately flush eyes or skin with a large amount of water for at least 15 minutes and seek immediate medical attention.

Quality Control:

The concentration of RNA should be determined by measuring the absorbance at 260 nm (A260) in a spectrophotometer. The ratio of the readings at 260 nm and 280 nm (A260/A280) provides an estimate of the purity of RNA. By RNX - Plus solution the ratio is usually greater than 1.6.

Signs

| Signs | Definitions | Signs | Definitions |
|---|----------------------------------|---|---|
|  | Temperature range on product use |  | Name and address of the manufacturer of the product |
|  | For Research Use Only |  | Product technical code |

Troubleshooting

This guide may help solve problems that may arise.

| Problem | Cause | Solution |
|-----------------------------|--|--|
| Lower than Expected 260/280 | Over dried pellet | Do not use speed vac. Dry pellet briefly at room temperature. |
| | Contamination of aqueous layer with interphase /organic phase | Take less aqueous phase. Use small bore pipette tips. Exercise care while removing aqueous layer. |
| | Sample contains glycogen, polysaccharides or other contaminants | Wash pellet in 4M LiCl prior to ethanol wash. |
| | Sample not homogenized with sufficient RNX Reagent. | Use 1 ml RNX Reagent for up to 50 mg tissue or 10 ⁶ cells. |
| Degraded RNA | Endogenous RNase Activity Exogenous | Use fresh tissue or cells. Process tissue immediately after removal. For cell culture samples, minimize washing steps. Add RNX Reagent directly to plates. |
| | RNase contamination | Store isolated RNA at -70°C, not -20°C. |
| | Improper storage of RNA | Extract samples within 20 minutes for multiple samples freeze homogenates at -70°C for later simultaneous processing. |
| | Homogenization Step extended beyond 20 minutes. | |
| DNA Contamination | Contamination of aqueous phase with interphase/organic phase | Take less of aqueous phase. Use small bore pipet tips. Exercise care while removing aqueous phase. |
| | Insoluble materials were not removed before chloroform extraction. | Remove any particulate material before chloroform addition. This material may trap DNA. |




شکایات مشتری

SinaClon
شرکت سیناکلون BioScience



فهرستی از مشتری

 Unit 9, Rouyesh building, Science and Technology Park, Tarbiat Modares University, Pajouhesh Blvd, Tehran, Iran



+982191082111



www.sinaclon.com



hi@sinaclon.com