

## dNTP Mix (Molecular biology grade)

**Concentration: 10mM**

**Cat. No.: MM2081 (Vol.:100μl)**

**Store at: -20°C**

**Cat. No.: MM2082 (Vol.:500 μl)**

**Shipment: Dry or Wet ice**

**Description:** Ready to use dNTP mixture. In this product concentration of each dNTPs is carefully balanced and formulated for higher and maximal fidelity of DNA synthesis.

**Composition:** SinaClon dNTP mix is prepared from individual ultrapure nucleotides. 10mM of each ultrapure dATP, dCTP, dGTP, dTTP in one vial and in aqueous solution (pH 7.0-7.5).

**Application:** For use in PCR, long PCR, RT-PCR, cDNA synthesis, primer extension, DNA sequencing and DNA labeling.

### Quality control

**Functional assay:** SinaClon dNTP mix was tested for cDNA synthesis, amplification of 977 and 788 bp multiplex PCR from human genomic DNA and DNA viruses (with SinaClon Taq DNA polymerase and Pfu polymerase).

### Individual Ultrapure Nucleotides

Molecular weight

1	dATP	557.2	C10H13N5O12P3Na3
2	dCTP	533.1	C9H13N3O13P3Na3
3	dGTP	573.2	C10H13N5O13P3 Na3
4	dTTP	548.1	C10H14N2O14P3Na3

### Basic PCR protocol

The following basic serves as a general guideline and a starting point for any PCR amplification. Optimal reaction conditions (incubation times, temperatures, concentration of Taq DNA polymerase, primer, MgCl<sub>2</sub> and template DNA) vary and need to be evaluated by the user. Add the following components to a sterile 0.5ml micro centrifuge tube sitting on ice:

Components	Volume	Final Concentration
<b>10XPCR Buffer (AMS)</b>	10μl	1X
<b>10mM dNTP mix</b>	2μl	0.2mM each
<b>50mM MgCl<sub>2</sub></b>	3μl	1.5mM
<b>Primers (10μM each)</b>	5μl	0.5μM each
<b>Template DNA</b>	1μg	-----
<b>Taq DNA polymerase</b>	0.5μl	2.5unit/100μl
<b>Autoclaved distilled water</b>	up to 100μl	-----

### PCR may be performed in 25-35 cycles as follows

<b>Denaturation</b>	<b>93°C</b>	45 seconds
<b>Annealing</b>	<b>55°C</b>	30 seconds
<b>Extension</b>	<b>72°C</b>	90 seconds

Optimal reaction conditions vary and need to be evaluated by the user. Mix and centrifuge dNTP before opening.



Unit 9, Rouyesh building, Science and Technology Park,  
Tarbiat Modares University, Pajouhesh Blvd, Tehran, Iran



**+982191082111**



**hi@sinaclon.com**



**www.sinaclon.com**

(NTP-056-00/00)