

SinaPure™ DNA

Kit for the isolation of DNA from Cell Culture, Tissue, Whole Blood & Gram Negative Bacteria

EX6011: 50 Preps

Store kit content at: RT

Keep enzymes at: -20°C

Keep tubes and bottles tightly closed.

Kit Contents	Quantity/ Volume	Unit	Storage Condition
Spin columns	50	PCS	RT
Collection tubes	50	PCS	RT
PreLysis solution	10	ml	RT
Proteinase K	1	ml	-20
Ribonuclease A	250	μl	-20
Lysis solution	10	ml	RT
Precipitation solution	10	ml	RT
Wash I solution	30	ml	RT
Wash II solution	50	ml	RT
Elution Buffer	5	ml	RT

Storage and Stability

The performance of the kit will be guaranteed until the printed expiration date. When the contents keep in the appropriate condition.

Kit Description

This kit contains the most necessary reagents for quick and Pure DNA preparation by spin column method from the recommended starting material:

Cell culture, Tissue, Gram Negative Bacteria and Whole Blood

SinaPure™ DNA kit system is one of the latest nucleic acid purification technologies. This kit presents remarkable features of timesaving, easy, prompt and high yield DNA purification. Basis of the technology is the binding of released DNA to matrices including silica membrane filters in presences of high salts concentration and reversibly elution in low salt condition. Obtained DNA is suitable for downstream applications including PCR.

Principle:

Samples would be degraded during incubation at the appropriate temperatures and in presence of PreLysis and Lysis solution. Meanwhile, enzymes accelerate the reaction by eliminating or reducing unwanted proteins and RNA contents of the cells simultaneously.

After successive washing steps, pure DNA would be released in low salt conditions.

Required Materials and equipment that were not provided:

- Heater block or water bath (37-56°C)
- Bench top micro centrifuge (12.000 RPM)
- Precision pipettes and sterile pipette tips (10 to 100 and up to 1000 μl).
- Sterile 1.5 ml or 2 ml polypropylene tubes

Important notes: please read before starting

- Follow GLP rules when performing nucleic acid extraction from biological samples.
- Check the procedure with a known sample, if use the kit for the first time.



- Check the vials for any possible crystal formation. In this case, place the vials at 37°C for 15 minutes and softly shake before use.
- Avoid freezing and thawing the samples.
- Use fresh samples.
- Keep the remaining initial samples preferably in small portions (200 µl, 20 mg...) or recommended volumes at minus -20 or colder.
- Aliquoting the enzymes may prolong their efficacy.

- Open spin columns only directly before use.

WARNINGS AND PRECAUTIONS

Caution: Avoid contact any kit reagents with skin & eyes. Wear PPE¹ before use any extraction kit. When handling biological samples, follow recommended procedures for biohazardous materials.

All steps should be performed at Room temperature.

Samples Preparation

Cell Culture

1. Collect cells ($1-3 \times 10^6$) by centrifugation at 5000 RPM for 5 minutes.
2. Dislodge the cell pellet by soft shaking or by finger tapping until it forms a complete suspension.
3. Add 200 µl PreLysis solution to the cell suspension.
4. Add 5 µl Ribonuclease A.
5. Vortex vigorously for 10 seconds.
6. Incubate at 37 °C for 15 minutes.
7. Add 20 µl Proteinase K.
8. Vortex vigorously for 10 seconds.
9. Incubate at 56°C for 30 minutes or until complete digestion.
Occasional (every 5-10 min.) and short vortex during incubation may enhance the final yield. In case of any remaining particles, centrifuge the tube briefly and transfer the upper phase to the new tube.
10. Add 200 µl Lysis solution to the tube and vortex vigorously for 10 seconds.
11. Incubate the tube at 56°C for 10 minutes.

Follow the laboratory protocol

Gram Negative Bacteria

1. Collect bacterial cells (up to 2×10^9) by centrifugation at 5000 RPM for 5 minutes.
2. Dislodge the cell pellet by soft shaking or by finger tapping until it forms a complete suspension.
3. Add 200 µl PreLysis solution to the cell suspension.
4. Add 5 µl Ribonuclease A.
5. Vortex vigorously for 10 seconds.
6. Incubate at 37 °C for 15 minutes.
7. Add 20 µl Proteinase K.
8. Vortex vigorously for 10 seconds.
9. Incubate at 56°C for 30 minutes or until complete digestion.
Occasional (every 5-10 min.) and short vortex during incubation may enhance the final yield. In case of any remaining particles, centrifuge the tube briefly and transfer the upper phase to the new tube.
10. Add 200 µl Lysis solution to the tube and vortex vigorously for 10 seconds.
11. Incubate the tube at 56°C for 10 minutes.

Follow the laboratory protocol

Blood

For the isolation of DNA from Whole blood, SinaPure™ Blood DNA, Cat.No. EX6001 is recommended.

Whole blood should be collected in EDTA (1mg/ml) to prevent blood clotting. Heparinized blood is not recommended.

Typically, 200µl of fresh blood is used for DNA isolation.

The blood can be kept at 2-8°C for (no longer than) 2 weeks for later use.

For long-term storage, samples should be aliquoted in 200µl portions and kept at -20°C.

To avoid any nuclease activity, keep samples frozen until DNA extraction.

1. Add 200 µl Lysis solution to the empty tube.
2. Add 20 µl Proteinase K and (optional) 5 µl Ribonuclease A.
3. Add 200 µl fresh blood to the tube.
4. Vortex vigorously for 20 seconds.
5. Incubate at 56°C for 15-30 minutes.

Follow the laboratory protocol

Tissue

Cut 5- 20 mg (maximum) of selected tissue and cut it into small pieces with a scalpel or homogenizer. Alternately, it can grind by mortar and pestle in liquid nitrogen. Freeze and thawed samples are not recommended. **(To cut the required amount of tissue refer to Tips and Suggestions Section).**

1. Add 200 µl PreLysis solution to the sample.
2. Add 20 µl Proteinase K and 5 µl Ribonuclease A.
3. Vortex vigorously for 10 seconds.
4. Incubate at 56°C for 1 to 3 hours (up to complete digestion).
 - For some tissues, incubation may continue overnight. In this case, complete digestion may need additional enzymes mix (20 µl).
 - Occasional (every 5-10 min.) and short vortex during incubation may enhance the final yield.
5. Add 200 µl Lysis solution to the tube and vortex vigorously for 10 seconds.
6. Incubate the tube at 56°C for 10 minutes.

In case of any remaining particles, centrifuge the tube briefly and transfer the upper phase to a new tube.

Follow the laboratory protocol

Protocol:

1	Add 200 µl Precipitation solution. Invert three times and transfer all tube content to the column quickly.
2	Centrifuge at 12,000 RPM for 1 min and Discard the filtrate.
3	Add 300 µl Wash I solution.
4	Centrifuge at 12,000 RPM for 1 min and discard the filtrate.



5	Add 300 µl wash I solution.
6	Centrifuge at 12,000 RPM for 1 min and discard the filtrate.
7	Add 700 µl wash II solution, Centrifuge at 12,000 RPM for 1 min and discard the filtrate.
8	Centrifuge empty column at 12,000 RPM for 2 min and discard the filtrate.
9	Place the tubes with the open lid into the heater block for 3- 5 min to evaporate any remaining ethanol.
10	Add 70 µl Elution Buffer directly at the center of the membrane and store at room temperature for 2-3 min, Centrifuge at 12,000 RPM for 1 min.

The final yield can be increased if:

Repeat centrifugation. Transfer eluted DNA to the center of the membrane once again and after a short incubation at room temperature (about 2 minutes) centrifuge at 12,000 RPM for 1 min.

Extend the incubation time in step 10 to about 10 min. (with closed lids).

You can discard the collection tube and place the column in a new 1.5 ml tube (not included).

Avoid contaminating the column with ethanol. Ensure that the column is dry and no ethanol contaminates the tip of the column. If you observe residues of ethanol, place the column in 56°C for 3-5 min.

Tips and Suggestions

Final elution volume depends on the initial sample: If higher DNA amounts are needed, a higher elution volume may use. Accordingly, a higher concentration may be achieved by a lower amount of Elution buffer.

Generally, 50-100 µl elution volume gives satisfactory results. An alternative way of increasing the DNA yield is repeated centrifugation. Transfer eluted DNA to the membrane filter once again and after a short incubation at room temperature (about 2 minutes) centrifuge at 12,000 RPM for 1 min.

For more purity you can add 500 µl wash II solution to the column in step 7 and centrifuge at 12,000 RPM for 1 min and repeat again. (provided)

Tissues: Final yields of DNA may vary depends on the different tissue types. In general, we found satisfactory results from 5 mg to maximum 20 mg of mice tissues like Kidney, Liver, Lung, Heart and Tail. However, in some cases you may need to optimize the initial starting materials (weigh 5- 20 mg) due to different levels of nucleic acids in specific tissue. Please consider, the final elution volume also needs to optimize for the appropriate concentration (50 – 150 µl).

Yields of genomic DNA will vary from sample to sample depending on the amount and type of material processed. In addition, the quality of starting material will affect DNA yield.

Troubleshooting

This guide may help solve problems that may arise.

Observation	Possible cause	Comments/suggestions
Low or no DNA yield	Inefficient lysis of sample	Make sure that: Incubation time with enzymes & homogenization step by vortex were enough. Decrease starting materials, extend incubation time and homogenization step. sample was frozen and thawed several times. Blood Sample was stored at 2-8°C longer than 2 weeks. Blood clots were present in the sample.
	Sample was frozen and thawed several times.	Keep samples freeze until DNA extraction. Whenever possible, fresh samples should be used and processed immediately. Storage at -20 °C is possible for several months. Several freeze-thaw cycles should be avoided, because this can result decreased molecular size of the DNA. Take new sample.
	Filter may clogged during purification	Check Lysis solution for any crystal formation. Warm Lysis before purification Check lysate for any tissues or particle remaining. Extend enzymatic incubation time, homogenization step or remove particles. Reduce the amount of initial sample.
DNA "smear"	Nuclease activity	Upon disintegration of Samples, cellular nucleases are released and may degrade genomic DNA. Use only sterilized glass and plastic ware in order to avoid nuclease Contamination.
Low DNA quality	Salt in elute	Repeat washing steps (I and II)
No enzymatic reaction	Residues of ethanol	Before adding Elution Buffer ensure that the column is dry, and no ethanol contaminates the tip of the column. If you observe residues of ethanol, place the column in a new reaction tube. Centrifuge again at 12,000 rpm for 1 min more

DNA quality control

- Agarose gel electrophoresis of prepared DNA is a direct method for testing DNA quality in terms of molecular size and conformation. Depending on the expected amount, pipette 3-5 ul eluted DNA directly to a gel slot. As for whole blood, DNA yield depends on quantity and quality of cells and storage duration and condition of sample.
- Photometric determination of DNA concentration and quality:
Determination of DNA concentration is done by UV reading at 260 nm. DNA preparations should be vortexed shortly and diluted (if needed) accordingly by 10 mM Tris-HCl or Elution buffer. Blank and dilution buffer should be the same. A standard procedure of

measuring DNA quality is the determination of the absorption quotient (Q) of readings at A260nm and A280nm:

$$Q = A260nm / A280nm$$

For a pure DNA preparation, Q lies between 1.7 and 2.0.

Kit Quality Control:

All components of this Kit are successfully tested in the DNA purification and amplification reaction for:
Tissue, cell Culture, gram negative bacteria and freeze or fresh whole blood.